

Quick Start Guide

The SpecPlate is a novel multiwell plate for plate-based absorption measurements and is designed to be integrated seamlessly into existing processes.

Nevertheless, this guide will help speed up its implementation at your lab.

Introduction The SpecPlate consists of 96 structures, each with an **inlet**, an **outlet**, and **four chambers** of descending height (2000 μm , 1400 μm , 700 μm , and 100 μm) for sample measurement (see figure 1). The **inlets follow a pattern of 96 wells**, and the samples flow clockwise from the inlet to the lowest chamber. The **air is displaced from the structures through the outlet**. The **chambers are arranged in a format of 384 wells** and can be measured using standard plate readers. A customized geometry is recommended, but measurements can also be conducted using pre-existing 384-well geometries. The **working volume per structure is 36 μL**

Compatibility The SpecPlate is **compatible with all common devices** capable of handling standard multiwell plates. This includes **liquid handlers** capable of handling 96-well plates and **plate readers** compatible with 384-well plates.

Manual Handling For best results, delicately place the pipette tip into the inlet and then **gradually and cautiously release** the solution into the structure. Do not dispense air gaps.

Automated handling The SpecPlate is designed to be compatible with well-calibrated liquid handlers, accommodating both **single and multi-channel pipettes** (8-channel, 96-channel). Please select appropriate pipette tips and flow rates for filling. Successful testing has been conducted with **flow rates ranging from 40 to 100 $\mu\text{L/s}$** , with a **pipetting height of about 2-3 mm** below the upper edge of the inlet opening. Slow lowering of the tips into the inlet is recommended. Handling with a **robotic manipulator arm** can be performed as usual. For stackers or similar devices, **adjust the necessary settings** in the labware definition file.

Disposal Disposal of used SpecPlates as usual with other plates.




Fig. 1 3D model of SpecPlate-structure




Fig. 2 Close-up of SpecPlate structure

Data Evaluation After measuring the SpecPlate in a 384-well format, for quick and easy data evaluation, it is recommended to **reorganize the measurement data by following the circular arrangement** of the measurement chambers associated with a measurement structure. This is illustrated schematically in the figure below. A simple Excel tool for this task can be provided on request.




Fig. 3 Scheme of data restructuring

Disclaimer Please note that PHABIOCC is not liable for damage to laboratory equipment caused by incorrect operation. SpecPlates are intended for **research use only**

Contact For further information visit www.phabioc.com or contact info@phabioc.com